

CLC2000™

For Needleless Access to Central, Peripheral, and Arterial Lines

The only IV connector guaranteed to eliminate the backflow of blood into the catheter lumen.



Female Luer



Male Spin Lock



PRIME the CLC2000

1. Remove from package and attach IV administration set or priming syringe to female luer of CLC2000.
2. Prime CLC2000 per facility protocol taking care to **invert and tap CLC2000 to expel air**.
3. Remove protective cap and attach male spin lock of CLC2000 to IV extension set or catheter.

SWAB and ACCESS the CLC2000

1. Swab female luer of CLC2000 before each access with desired disinfectant per facility protocol.
2. To access with luer lock or luer slip IV device, ensure that the **luer is fully prime** to prevent the introduction of air.
3. Grasping CLC2000, firmly push and twist connecting device into CLC2000 until tight. To remove device, grasp CLC2000 and twist device away until loose.

NOTE: If using a spin lock device, first secure the luer and then "lock down" spin collar. This will ensure a secure connection and optimal flow rates.



FLUSH the CLC2000

1. Flush the CLC2000 after each use per facility protocol.
2. **DO NOT** clamp tubing or catheter while flushing or disconnecting. Clamping may interrupt positive displacement.
3. The CLC2000 is a closed barrier to bacteria. **Do not use end caps.**

CAUTION:

- Federal (USA) law restricts this device to sale by or on the order of a physician.
- It is recommended this device be changed per CDC guidelines or per validated facility protocol.

