

# Not **All** Feeding Tubes Are Created Equal!



# CORFLO<sup>®</sup>

**CORFLO<sup>®</sup>** Enteral Feeding Tubes from VIASYS MedSystems offer advantages to your patients **unlike any other feeding tubes available**, including:

● Patented **ULTRA** anti-clog outlet port significantly reduces the occurrence of clogging.

● The largest selection of adult, pediatric, and neonatal feeding tubes to allow one-stop sourcing for all your feeding tubes hospital or system-wide.

● Patented **TRAK** feeding tubes: The **ONLY** feeding tubes that can be used with the **CORTRAK<sup>™</sup>** Enteral Access System.

● **CORFLO Anti-IV** prevents inadvertent feeding into an IV line. Compatible only with oral tip syringes for neonatal and pediatric use.



# The CORFLO<sup>®</sup> Advantage

## The Ultimate in Patient Safety:

- All CORFLO feeding tubes are printed with centimeter marks for increased patient safety and intubation accuracy.
- The entire CORFLO polyurethane feeding tube line is color coded to correspond to varying lengths.
- Stylet lengths are matched to the CORFLO feeding tubes.

## Clinician Convenience and Patient Comfort:

- Flow-through stylets to flush, auscultate and aspirate during the insertion procedure
- Water activated lubricant on all polyurethane feeding tubes for greater patient comfort during intubation and effortless stylet removal
- 100% radiopaque for X-ray confirmation
- Weighted tubes use multiple tungsten weights allowing for flexibility and greater comfort during intubation
- "Y" ports on all CORFLO polyurethane feeding tubes minimize touch contamination and makes flushing easy

*VIASYS MedSystems feeding tubes are the standard by which all other feeding tubes are measured!*

VIASYS MEDSYSTEMS 20-9431 10 Fr 43"-109cm 100 + + + + 95 +

FR Size/ Length	5 FR STERILE	6 FR	8 FR	10 FR	12 FR	
15" (38cm)	20-1155AIV* 20-1155* 20-PVC155*	20-1156AIV* 20-1156*	20-1158AIV* 20-PVC158			<b>CORFLO<sup>®</sup> Anti-IV:</b> developed to prevent tubing misconnections for pediatric & neonatal use; will only accept ORAL Syringes. Sterile & DEHP FREE*** <b>CORFLO-ULTRA NG**:</b> designed for intubations where anti-clogging and ease of insertion are critical. <b>CORFLO-CONTROLLER NG:</b> developed for maximum control during the intubation procedure especially with Endoscopy and Fluoroscopy placements. <b>CORFLO-INFANT PVC NG:</b> manufactured from premium polyvinyl chloride, packaged sterile for short-term feeding for infants. <b>CORFLO-ULTRA Lite NG</b> Non-wt. w/stylet <b>CORFLO-ULTRA Lite NG</b> Non-wt. w/o stylet <b>CORFLO-ULTRA Lite Clear NG</b> STERILE, Clear, non-wt. w/o stylet <b>CORFLO-ULTRA Pedi NG</b> 1 g & 1.5g Pediatric wt. w/stylet <b>CORFLO-ULTRA NG</b> 3 g wt. w/stylet <b>CORFLO-ULTRA NG</b> 3 g wt. w/o stylet <b>CORFLO-ULTRA Pili NG</b> 3 g wt. w/stylet <b>CORFLO-ULTRA 7 NG</b> 7 g wt. w/stylet <b>CORFLO-CONTROLLER NG</b> 3 g wt. w/stylet <b>CORFLO-CONTROLLER Pili NG</b> 3 g wt. w/stylet <b>CORFLO-CONTROLLER 7 NG</b> 7 g wt. w/stylet <b>CORFLO-INFANT P NG</b> STERILE, Plastic, non-wt. w/o stylet <b>CORTRAK ENTERAL ACCESS SYSTEM NG</b> Non-wt. w/stylet <b>CORFLO Anti-IV NG</b> STERILE, Clear, non-wt. w/o stylet *STERILE **All CORFLO-Ultra feeding tubes have our unique ANTI-CLOG feeding port. ***Contact VIASYS MedSystems for Extension Set re-order numbers or for further information regarding CORFLO Anti-IV. Boluses: approximate weights
22" (56cm)	20-1225AIV* 20-9225* 20-C1225*	20-9226 20-1226 20-1226AIV* 20-8226 20-C1226	20-1228AIV* 20-9228 20-C1228			
36" (91cm)	20-1365*	20-9366 20-1366 20-8366 20-7366	20-9368 20-1368 20-8368 20-7368	20-9361 20-1361 20-7361 20-5361	20-1362 20-5362	
43" (109cm)			20-9438TRAK 20-9438 20-1438 20-7438 20-5438 20-6438 20-4438 20-3438 20-2438 20-0438	20-9431TRAK 20-9431 20-1431 20-7431 20-5431 20-6431 20-4431 20-3431 20-0431 20-2431	20-9432 20-1432 20-7432 20-5432 20-4432 20-0432	
55" (140cm)	All Feeding Tubes packaged 10/case, except Infant PVC (50/case)		20-9558	20-9551 20-9551TRAK 20-3551 20-2551		

**Enterally speaking, nobody delivers better!®**

VIASYS MedSystems  
100 Chaddick Dr.  
Wheeling, IL 60090

1-800-323-6305  
1-847-403-3400 • FAX 1-847-541-9526  
www.viasyshealthcare.com

VIASYS  
HEALTHCARE  
MEDSYSTEMS

CE 0144  
ISO 13485  
5811E-09-06