

CORFLO®

Anti-IV

Enteral Feeding System

**A new enteral feeding standard of care
for uncompromised patient safety**

For Neonatal and Pediatric Patients



- **Unique design is compatible ONLY with oral tip syringes or fittings**
- **Proprietary feeding port with CORFLO Anti-IV male connector**



CORFLO Anti-IV is a complete enteral feeding system that:

- Prevents inadvertent feeding into an intravenous line
- Prevents injecting intravenous medication into the feeding port of the CORFLO Anti-IV feeding tube
- Prevents direct connection of the feeding set or NG tube to any standard luer lock or luer slip syringe
- Fits all standard oral tip syringes and adapted syringes with oral tips
- Meets JCAHO Sentinel Event Alert Recommendations (Issue 36, April 3, 2006)



CORFLO® *Anti-IV* System Product Specifications

CORFLO Anti IV Feeding Tubes

	Length	French	#/case
20-1155AIV	15"	5FR	10
20-1156AIV	15"	6FR	10
20-1158AIV	15"	8FR	10
20-1225AIV	22"	5FR	10
20-1226AIV	22"	6FR	10
20-1228AIV	22"	8FR	10
20-1365AIV	36"	5FR	10
20-1366AIV	36"	6FR	10
20-1368AIV	36"	8FR	10

CORFLO Anti-IV Extension Sets

	Description	Approx. Priming Volume	#/case
20-1015AIV	15" Minibore	0.5 ml	50
20-1030AIV	30" Minibore	1.1 ml	50
20-1060AIV	60" Minibore	2.25 ml	50
20-1060AIV-S	60" Minibore w/ Step Connector	2.25 ml	50
20-1030AIV-LB-S	30" Large Bore w/ Step Connector	4.0 ml	50
20-1006AIV-LB	6" Large Bore	0.8 ml	50
20-1030AIV-LB	30" Large Bore	4.0 ml	50

The CORFLO Anti-IV Feeding Tube and Extension Set can be used with any feeding method. The proprietary connection on the feeding tube will allow feeding only through the CORFLO Anti-IV Feeding Set or ORAL tip syringe to insure enteral only use!

FEEDING TUBE:

- Medical grade polyurethane
- Color: Clear
- Radiopaque orange stripe
- Product identification printed on tube
- Centimeter markings (approx.) printed on tube to aid placement and check migration
- Unique "Y" Port—No disconnecting required and reduces touch contamination
- Feeding Port: Proprietary CORFLO Anti-IV Male Luer Connector
- Sterile
- DEHP Free
- Latex Free

FEEDING EXTENSION SETS:

- Compatible with any size oral syringe
- Fittings to prevent inadvertent administration into IV, PICC or Central Line
- Sterile
- DEHP Free
- Latex Free



Oral Syringe Ordering Options

Ameritus Medical Products: 1-800-825-5996

www.kentecmedical.com

Baxa Corporation: 1-800-395-2292

www.baxa.com

Becton Dickinson (BD): 1-888-237-2762

www.bd.com

Monoject: 1-800-962-9888

www.kendall.com

CORPAK MedSystems

100 Chaddick Dr., Wheeling, IL 60090

Call Toll Free: 1-800-323-6305

Phone: 1-847-403-3400 Fax: 1-847-541-9526

www.corpakmedsystems.com

www.corfloanti-iv.com