

CORFLO[®] Feeding Tube Matrix

CORFLO[®] Anti-IV: developed to prevent tubing misconnections for pediatric & neonatal use; will only accept ORAL Syringes. Sterile & DEHP FREE***

CORFLO-ULTRA NG:** designed for intubations where anti-clogging and ease of insertion are critical.

CORFLO-CONTROLLER NG: developed for maximum control during the intubation procedure especially with Endoscopic and Fluoroscopic placements.

FR Size/ Length	5 FR Sterile	6 FR	8 FR	10 FR	12 FR
15" (38 cm)	20-1155AIV*	20-1156AIV*	20-1158AIV*		
	20-1155*	20-1156*			
22" (56 cm)	20-1225AIV*	20-9226	20-1228AIV*		
	20-9225*	20-1226	20-9228		
	20-C1225*	20-1226AIV*	20-C1228		
		20-8226			
		20-C1226			
36" (91 cm)	20-1365*	20-9366	20-9368TRAK	20-9361TRAK	20-1362
	20-1365AIV*	20-1366	20-9368	20-9361	20-5362
		20-8366	20-1368	20-1361	
		20-7366	20-8368	20-7361	
		20-1366AIV*	20-7368	20-5361	
			20-1368AIV*		
43" (109 cm)			20-9438TRAK	20-9431TRAK	20-9432
			20-9438	20-9431	20-1432
			20-9438AIV	20-9431AIV	20-7432
			20-1438	20-1431	20-5432
			20-7438	20-7431	20-4432
			20-7438AIV	20-7431AIV	20-0432
			20-5438	20-5431	
			20-6438	20-6431	
			20-4438	20-4431	
			20-3438	20-4431AIV	
			20-2438	20-3431	
			20-2438AIV	20-0431	
			20-0438	20-0431AIV	
				20-2431	
55" (140 cm)			20-9558	20-9551	
				20-9551TRAK	
All Feeding Tubes packaged 10/case.				20-3551	
				20-2551	
1.5 g wt.					

CORFLO-ULTRA Lite NG	Non-wt. w/stylet
CORFLO-ULTRA Lite NG	Non-wt. w/o stylet
CORFLO-ULTRA Lite Clear NG	STERILE, Clear, Non-wt. w/o stylet
CORFLO-ULTRA Pedi NG	1 g & 1.5 g Pediatric wt. w/stylet
CORFLO-ULTRA NG	3 g wt. w/stylet
CORFLO-ULTRA NG	3 g wt. w/o stylet
CORFLO-ULTRA Pili NG	3 g wt. w/o stylet
CORFLO-ULTRA 7 NG	7 g wt. w/stylet
CORFLO-CONTROLLER NG	3 g wt. w/stylet
CORFLO-CONTROLLER Pili NG	3 g wt. w/stylet
CORFLO-CONTROLLER 7 NG	7 g wt. w/stylet
CORTRAK ENTERAL ACCESS SYSTEM NG	Non-wt. w/stylet
CORFLO Anti-IV NG	STERILE, Clear, non-wt. w/o stylet

*STERILE

**All CORFLO-Ultra feeding tubes have our unique ANTI-CLOG feeding port.

***Contact CORPAK MedSystems for CORFLO Anti-IV Extension Set re-order numbers or for further information regarding CORFLO Anti-IV.

Boluses: approximate weights