





Nutri-Safe and Solu-Safe Extension Sets – Groups C and G

Catalog Number	Description	Length (in-cm)	Priming Volume	Unit per Case
27117	24" NutriSafe™-C Extension Set with 1.2 μ filter, one CLAVE® Needle-Free "Y" port and spin luer lock. 	24" - 60cm	5 ml	150
27118	16" Solu-Safe™-C Extension Set with 0.22 μ filter, one CLAVE® Needle-Free "Y" port, and spin luer lock. 	16" - 40cm	4 ml	150
27074	24" Nutri-Safe™-G Extension Set with 1.2 μ filter, one latex-free standard "Y" port and spin luer lock. 	24" - 60cm	5 ml	150
27075	16" Solu-Safe™-G Extension Set with 0.22 μ filter, one latex-free standard "Y" port and spin luer lock. 	16" - 40cm	4 ml	150