







Victus Universal Standardbore IV Sets – Group C with CLAVE® Needle Free Ports

Catalog Number	Description	Length (in-cm)	Priming Volume	Unit per Case
27071	74" BasIVset™-C Vented I.V. Set with macrodropper (20 drops/ml), one CLAVE® Needle-Free "Y" port and spin luer lock. 	74" - 187cm	15 ml	100
27114	82" Microset™-C Vented I.V. Set with microdropper (60 drops/ml), one CLAVE® Needle-Free "Y" port, and spin luer lock. 	82" - 208cm	15 ml	100
27072	82" Microset™-C Vented I.V. Set with macrodropper (20 drops/ml), two CLAVE® Needle-Free "Y" ports and spin luer lock. 	82" - 208cm	15 ml	100
27080	87" UniVset™-C Vented I.V. Set with macrodropper (20 drops/ml), two CLAVE® Needle-Free "Y" ports, check valve and spin luer lock. 	87" - 220cm	16 ml	100
27115	92" Bureset™-C I.V. Set with calibrated burette (100 ml), with CLAVE® injection port automatic shut-off floating valve, one CLAVE® Needle-Free "Y" port, and spin luer lock. 	92" - 234cm	15 ml	50
27116	92" Bureset™-C I.V. Set with calibrated burette (150 ml), with CLAVE® injection port automatic shut-off floating valve, one CLAVE® Needle-Free "Y" port, and spin luer lock. 	92" - 234cm	15 ml	50