







CLAVE Needle free system – Closed Vial Access Spikes for Chemotherapy

Catalog Number	Description	Image	Units per Case
CH-77	GENIE Multi-Dose vial spike with pressure equalizing internal balloon system. Stability skirt to fit any 20 mm vial, access through CLAVE closed connector		50 per case
CH-70	Dual Vented Universal Multi-Dose vial spike. Lateral fittings to attach to all size vials, access through CLAVE closed connector		50 per case
CH-62	Multi-Dose vial spike. Stability skirt to fit any 14 mm vials, access through CLAVE closed connector		50 per case
CH-50	Non Vented Universal Multi-Dose vial spike. No support, attaches to all size vials, access through CLAVE closed connector		50 per case
CH-71	Dual Vented Universal Multi-Dose vial access spike. Lateral fittings to attach to all size vials.		50 per case
CH-60	Universal single dose vial access spike, non vented, no support.		50 per case