







CLAVE Needle free system– Closed Chemotherapy Administration IV Sets

Catalog Number	Description	Units per Case
CH3065	15 drop administration set, rotating Luer connector, with hanger, 41" (104 cm), Non-DEHP tubing 	50 per case
CH3025	15 drop administration set, SPIROS connector, with hanger, 41" (104 cm), Non-DEHP tubing 	50 per case
CH3072	20 drop CLAVE chamber administration set, rotating Luer connector, with hanger, 30" (76 cm), Non-DEHP tubing 	50 per case
CH3011	20 drop CLAVE chamber administration set, SPIROS connector, with hanger, 30" (76 cm), Non-DEHP tubing 	50 per case
CH3018	20 drop administration set, with proximal SPIROS connector, locking cannula, 40" (102 cm), Non-DEHP tubing 	50 per case
CH3078	Specialty extension / administration set with two CLAVE connectors, spike and female spike adapter, 12" (30 cm), Non-DEHP tubing 	50 per case