








## Victus Universal Standardbore IV Sets – Group G

Catalog Number	Description	Length (in-cm)	Priming Volume	Unit per case
27058	74" BasIVset™-G Vented I.V. Set with macrodropper (20 drops/ml), one latex-free standard "Y" port and spin luer lock 	74" - 187cm	15 ml	100
27060	74" Microset™-G Vented I.V. Set with microdropper (60 drops/ml), one latex-free standard "Y" port, and spin luer lock 	74" - 187cm	15 ml	100
27059	81" Macroset™-G Vented I.V. Set with macrodropper (20 drops/ml), two latex-free standard "Y" ports and spin luer lock 	81" - 207cm	15 ml	100
27043	81" UniVset™-G Vented I.V. Set with macrodropper (20 drops/ml), two latex-free standard "Y" ports, check valve and spin luer lock 	81" - 207cm	16 ml	100
27062	92" Bureset™-G I.V. Set with calibrated burette (100 ml), automatic shut-off floating valve, one latex-free standard "Y" port, and spin luer lock 	92" - 234 cm	15 ml	50
27063	92" Bureset™-G I.V. Set with calibrated burette (150 ml), automatic shut-off floating valve, one latex-free standard "Y" port, and spin luer lock 	92" - 234cm	15 ml	50
27061	36" Secunset™-G Vented I.V. Set with macrodropper (20 drops/ml) and spin luer lock 	36" - 92cm	9 ml	100